



Available only at
trioSKINCARE
Anti Aging ▲ Acne Management Clinic

What is SkinMASTER™? An Ultrasonic Skin Rejuvenating System beyond microdermabrasion..

A new patented technology that provides non-invasive solutions to diminish fine lines, wrinkles, dehydration, blemishes, enlarged pores, irregular skin texture and skin discolorations. The SkinMASTER incorporates multiple technologies with six unique precision instruments for a variety of customized treatments to meet your specific skin needs. Say goodbye to microdermabrasion and hello ultrasonic skin rejuvenation.

How does the SkinMASTER™ work? In just 2 steps you can have vibrantly different skin!

1. Ultrasonic hydro-cavitation: A series of vibrations are transmitted to the skin and underlying tissue via elastic waves and water at frequencies as high as 25,000 cycles per second to create both mechanical and thermal action for an effective peeling process. A unique hand piece is used to ultrasonically exfoliate and deeply cleanse the outer layer of the skin. This is the first step in all SkinMASTER treatments to effectively remove dead skin cell build up and reveal vibrant healthy skin. Exfoliation intensity is controlled by your skin care specialist and will vary with each individual's skin type and condition.

2. Ultrasonic Elastic Wave: This is a low frequency ultrasound that micronizes or breaks down an ingredient into tiny particles so that your skin can absorb and utilize the vitamins and nutrients efficiently and effectively. The "micro massage" action of elastic waves dramatically increases the skin hydration levels while promoting deep relaxation and skin plumping. This deep penetrating moisture therapy is included in all SkinMASTER treatments. The therapeutic product used in this step will be customized to meet the individual's skin type and condition.

Options:

3. Modified Lotti Micro-Technology (MLM): A disc shaped wand delivers electrical stimulation to re-educate facial muscles and dramatically increase skin circulation. MLM is used to promote increased skin health by stimulating blood flow and increased oxygen levels in the skin. MLM lifts and firms facial contours by promoting muscular contractions and strengthening. Ask your skin care specialist if this technology would benefit your skin condition and skin health.

4. Therapeutic Microdiode: This is the delivery of amplified light energy not only to increase blood circulation but to stimulate the production of new healthy collagen and elastin. Microdiode therapy can immediately plump fine lines around the eyes, mouth, lips and forehead while regenerating collagen and elastin in the deep dermal areas of the skin. Slight discomfort, erythema and mild swelling are associated with microdiode treatment.

5. Ultrasonic Blemish Clearing: The SkinMASTER treatment can be intensified by utilizing a clearing module which combines ultrasonic cavitation and elastic wave therapy to aggressively exfoliate and then penetrate melanin regulating agents (bleaching agents) to accelerate the clearing of skin discolorations.

More Benefits to SkinMASTER™

- Helps diminish bacteria and excess oil to purify skin, reduce blemishes and reduce pore size.
- Leaves skin hydrated, oxygenated and rejuvenated.
- Excellent for even the most sensitive skin types – No chemicals, just ultrasound and water!
- No redness, No irritation, No crystal residue, No downtime.
- While all SkinMASTER treatments include extraction, your skin care specialist will advise you if an alternative treatment would be more effective for your skin condition and concerns.
- Customize your treatment to focus on your key areas of concern. Starting @ \$125.00

Traditional microdermabrasion is limited to Exfoliation only. With SkinMASTER you can Exfoliate, Hydrate, Penetrate, Circulate, Plump, and Tighten the skin while provide Lymphatic Drainage.

MICRODERMABRASION

The Past and Present

History

Dermabrasion therapies are not a recent discovery. Traces of abrasion techniques were discovered among the Egyptians as a means to obtain smoother skin with abrasive materials. In the early 1900's physicians utilized static abrasive disc combined with suction to polish the skin. This process has been refined through present day and is performed under anaesthesia. It causes painful bleeding, scabbing and irritation for a minimum of 4 weeks. From this technique the design of the particle microdermabrasion system emerged. The MicroDermabrasion system was developed by Mattioli Engineering of Florence, Italy. It has been utilized in Europe for many years prior to its FDA approval in the United States in 1996.

What is Particle Microdermabrasion?

The system uses an intense buffing process which resurfaces the skin by utilizing fine aluminum oxide crystals to remove the upper layer of the stratum corneum. At the same time the crystals are projected onto the skin, the hand piece of the instrument contains a vacuum tube to suction up both crystals and the skin's cellular debris.

What can it treat?

This process exfoliates the epidermis, resulting in refined skin texture and reduced pore size. The skin appears brighter, softer and smoother.

What can it *not* treat?

It cannot treat acne, hyperpigmentation, rosacea or wrinkles. It cannot treat dehydration.

Controversies -

1. Aluminum Oxide Crystals
2. Manufacturer Claims - Effective acne and pigmentation treatment
3. Exasperation of other skin conditions.

trioSkincare has gone beyond simple microdermabrasion and introduce Ultrasonic *particle-free* Skin Rejuvenation.

What is Ultrasonic Skin Rejuvenation?

A system that uniquely incorporates 5 technologies without the irritation, inflammation and risks associated with traditional particle abrasion machines. It is a system that provides benefits for acne, pigmentation, rosacea, oily, dry, aging and sensitive skin.

What can it Do?

1. Ultrasonic Hydro-Cavitation. - Exfoliation
2. Sonophoresis - Hydration
3. Amplified Light Energy – Class 111A laser that delivers amplified light energy to stimulate circulation in the dermis, increase collagen and elastin production.
4. Modified Lotti Micro-Technology (MLM) -electrical stimulation of facial muscles for increased collagen and elastin production.
5. Hyper-Oscillating Vibration – High frequency of hand piece creates a complete antibacterial environment.